


Dosierventil ST-161

Dosierventil ST-161 Ansaugverhalten					
Skalenwerte	1,3 - 1,7 (10 l/min)	1,6 - 2,0 (15 l/min)	1,8 - 2,3 (20 l/min)	2,4 - 2,8 (30 l/min)	2,8 - 3,2 (45 l/min)
1	0,139	0,106	0,220	0,198	0,210
2	0,430	0,446	0,480	0,250	0,460
3	0,490	0,550	0,582	0,600	0,584
4	0,618	0,688	0,744	0,836	0,722
5	0,656	0,740	0,876	1,012	0,856
6	0,680	0,802	0,924	1,090	0,902
7			1,030	1,148	1,012
0			1,090	1,220	1,060
1				1,374	
2					
3					
4					
5					
6					
7	0,726	0,916	1,160	1,608	1,112

Dosierventil ST-161



Dosierung mit Beschränkungen/Dosiergrößen

Durchmesser Dosiergröße	1,3 - 1,7 (10 l/min)	1,6 - 2,0 (15 l/min)	1,8 - 2,3 (20 l/min)	2,4 - 2,8 (30 l/min)	2,8 - 3,2 (45 l/min)
0,5	0,202	0,318	0,228	0,220	0,210
0,6	0,348	0,330	0,232	0,320	0,215
0,7	0,398	0,345	0,414	0,444	0,410
1,0	0,558	0,652	0,676	0,678	0,660
1,1	0,604	0,714	0,744	0,736	0,720
1,2	0,615	0,750	0,822	0,810	0,798
1,3	0,646	0,806	0,958	0,938	0,930
1,5	0,662	0,850	1,076	1,172	1,045
ohne Düse	0,732	0,942	1,260	1,566	1,190

